

**Butterfly Blood Collection needle**

- Simple safety
- Reduces the risk of
- Thumb rest to ease handling during safety activation
- Flexible wings to secure device to patient's skin
- Thin walled needle is designed to minimize patient discomfort
- Available in 21, 23 or 25 gauge
- Transparent housing and tubing for clear visualization of blood flow
- DEHP-free, flexible, soft, and kink resistant tubing



| Description                       | Cat. No.  | Color Code | Size G/" | Diameter | Length |
|-----------------------------------|-----------|------------|----------|----------|--------|
| Butterfly Blood Collection needle | CY1102725 |            | 27G x 1" | 0.40 mm  | 13 mm  |
|                                   | CY1102525 |            | 25G x 1" | 0.50 mm  | 19 mm  |
|                                   | CY1102425 |            | 24G x 1" | 0.55 mm  | 19 mm  |
|                                   | CY1102325 |            | 23G x 1" | 0.60 mm  | 25 mm  |
|                                   | CY1102225 |            | 22G x 1" | 0.70 mm  | 25 mm  |
|                                   | CY1102125 |            | 21G x 1" | 0.80 mm  | 27 mm  |
|                                   | CY1102025 |            | 20G x 1" | 0.90 mm  | 28 mm  |
|                                   | CY1101925 |            | 19G x 1" | 1.20 mm  | 28 mm  |

**\* Packing info**

| Product                           | Needle Size | Primary       | Secondary    | Packing Size | G.W.    | N.W.    |
|-----------------------------------|-------------|---------------|--------------|--------------|---------|---------|
| Butterfly Blood Collection needle | 20G to 27G  | 100 pcs / bag | 50 bag / Ctn | 65*40*37     | 13.6 Kg | 11.9 Kg |